



## READY LEMON DISINFECTANT

**Kills SARS-CoV-22 virus, the cause of COVID-19. Effective against SARS-Related Coronavirus 22 [SARS-CoV-2] in 3 minutes.**



Use **LEMON DISINFECTANT** as a disinfectant on hard, nonporous surfaces. Use for non-scratch cleaning of showers and tubs, shower doors and curtains, fixtures and toilet bowls. Use to clean and disinfect finished floors, sinks and tubs. Use with a mop and bucket, trigger sprayer, sponge or by soaking.

Use **LEMON DISINFECTANT** in: Health care facilities, Restaurants [Front of House], Schools, and Institutional facilities. Cross-contamination is of major housekeeping concern. **LEMON DISINFECTANT** has been formulated to aid in the reduction of cross-contamination not only in hospitals, but in schools, institutions and industry.

Use **LEMON DISINFECTANT** on washable hard, nonporous surfaces of: bathroom fixtures, floors, and non-critical medical equipment surfaces: medical examining tables, medical lamps, rescue tools, resuscitators, medical scales, stands, stretchers, stethoscopes, walkers and wheelchairs. Other hard nonporous surfaces made of: plastic and stainless steel.

### DIRECTIONS FOR USE

#### DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Before use in federally inspected meat and poultry food processing plants and dairies, food products and packaging materials must be removed from the room or carefully protected.

#### DISINFECTION / VIRUCIDAL\* / FUNGICIDAL / MOLD & MILDEW CONTROL DIRECTIONS:

Spray 6 - 8 inches from surface, until surfaces are thoroughly wet. Do not breathe spray. Treated surfaces must remain wet for 10 minutes. **For Influenza Virus Type A and Human Coronavirus**, treated surfaces must remain wet for 1 minute. **Effective against SARS-CoV-2, the virus that causes COVID-19, in 3 minutes.** Wipe dry with a clean cloth, sponge, or mop or allow to air dry. For heavily soiled areas, a preliminary cleaning is required.

Rinse all surfaces that come in contact with food such as countertops, appliances, tables and stovetops with potable water before reuse.

**For fungicidal activity: LEMON DISINFECTANT** is an effective one-step fungicide against *Trichophyton mentagrophytes* when used on hard nonporous surfaces in areas such as locker rooms, dressing rooms, shower and bath areas and exercise facilities. Allow surface to remain wet for 10 minutes.

**For mold and mildew:** Repeat treatment every 7 days or more often if new growth appears.

SEE PRODUCT LABELING FOR MORE USAGE DIRECTIONS FOR SPECIFIC APPLICATIONS.

*\*Please consult your ACS sales Representative or call 1-800-889-3489 for additional use directions.*

### GENERAL PRODUCT INFORMATION

**Description:** Ready to Use Disinfectant, Cleaner, Mildewstat, Fungicide, Deodorizer, Virucide

**Safety:** Keep product stored in original container, in an upright position, with container tightly closed. If swallowed, rinse mouth - Do Not induce vomiting. If on skin, wash with water - remove contaminated clothing. If inhaled, remove person to fresh air - administer artificial respiration if not breathing. If in eye, flush eye continuously with water - remove contacts and continue flushing with water. Seek medical assistance if symptoms persist. Dispose of contents/container according to local/regional regulations. For Chemical Emergency, spill, leak, fire, exposure, or accident, contact: Chemtrec (day or night) - Within USA and Canada - 1-800-424-9300.

**Packaging:** 2.75 Gal. (10.41L)

**Product #:** 19001903

See Safety Data Sheet (SDS) for a complete listing of hazards, precautions, first aid statements, storage and handling information.

### AMERICAN CHEMICAL SYSTEMS



Manufactured by:

**AMERICAN CHEMICAL SYSTEMS II, LLC**

3505 W. 29th Street South

Wichita, Kansas 67217 • 800-889-3489

<http://americanchemicalsystems.com/>

READY PRODUCTS